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Confidential - distribution prohibited without permission.

Consistency - a critical fact in IVD product development

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
Laura-Leena Kiiskinen, PhD
R&D Director, Bioreagents

Confidential - distribution prohibited without permission.

March 18 2018


Consistency 一致性

"The quality of achieving a level of performance which does not vary greatly in quality over time."

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-Definition of consistency in English by Oxford Dictionaries-
-牛津字典中一致性的英语定义-

Confidential - distribution prohibited without permission.


Antibody consistency 抗体的一致性

- Antibodies ultimately determine the performance of an IVD assay
抗体根本上决定了IVD检测的性能
- Each antibody has its unique specificity, affinity, surface charge, and solubility properties
每个抗体都各自具有唯一的特异性、亲和力、表面电荷及溶解度。
- Consistent quality means that the properties do not change from batch to batch
质量的一致性意味着不同批号产品的性能不会发生变化

Confidential - distribution prohibited without permission.

Critical antibody features 关键抗体性能

Confidential - distribution prohibited without permission.


Example: Consequences of antibody inconsistency 抗体一致性不佳的后果

Feature 特性	Consequences of inconsistency 抗体一致性不佳的后果				
	Re-optimization of antibody amount for each new IVD product lot 须根据每批IVD产品重新优化抗体用量	Sensitivity problems 灵敏度问题	Need to optimize the labelling/coating process pH or volume 须优化包被标记过程的PH及体积	Removal of aggregates, fragments, or impurities needed 须去除聚合物、片段或杂质	Higher antibody amounts needed in the process 过程中须使用更多的抗体
Immunoreactivity 免疫学活性	✓	✓			✓
Homogeneity 同质性	✓	✓		✓	✓
Charge isoform profile 电荷亚型分布		✓	✓		
Purity 纯度	✓	✓		✓	✓
Concentration 浓度		✓	✓		

Confidential - distribution prohibited without permission.

Antibody consistency plays a role in all stages of an IVD product' s life cycle

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- Product development 产品开发
- Product validations 产品确认
- Manufacturing 产品生产

Confidential - distribution prohibited without permission.

Antibody consistency in an IVD product' s life cycle IVD产品生命周期中的抗体一致性

IVD product life cycle phase IVD产品生命周期的阶段	Effects caused by poor antibody consistency 抗体一致性差的影响	Consequences 后果
Product development 产品开发	Inconclusive feasibility results, repeated testing 不确定是否可行的结果，反复测试	Extra time and cost, project rejection 额外的时间及费用，拒绝立项
Product validation 产品确认	Failure in reproducibility testing conducted with different raw material batches* 采用不同批次原材料进行重复测试时测试失败	Extra time and cost, project discontinuation, loss of sales 额外的时间及费用，项目取消，销售下降
Manufacturing 产品生产	Lot-to-lot variation in the final IVD product 最终产品批间差异	Scrap costs, extra QC costs, loss of sales 报废成本，额外质检费用，销售下降

Confidential - distribution prohibited without permission.

*common industry practice: ≥ 2 raw material batches must be used in the validation testing of a new IVD product

How does Medix Biochemica ensure antibody consistency?

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- ISO 13485 certified quality management system

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- Precise and comprehensive quality control analyses

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- Focus on IVD market requirements in all company functions, from product ideas to product deliveries

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Confidential - distribution prohibited without permission.

 Medix Biochemica

High antibody quality begins in R&D 高质量的抗体始于研发

Product idea generation 产品方案的产生

- Ideas collected from IVD customers
从IVD客户中收集创意
- Prioritization based on customer needs and technical feasibility
根据客户需求及技术可行性建立研发优先级

Product development 产品开发

- Patient samples tested early in the development process
研发初期就进行患者样本测试
- Thorough biochemical and immunological characterization of new mAbs对新抗体进行全面的生化及免疫学
- The best cell lines adapted to serum-free culture medium, to ensure reproducibility
采用最好的细胞株进行无血清培养，从而保证可重复性

Product validation 抗体确认

- Antibody production, purification, and quality control performed in a similar manner as in manufacturing 采用与生产相似的抗体生产、纯化、质量控制方法

Confidential - distribution prohibited without permission.

 Medix Biochemica

Analyses performed in R&D for new monoclonal antibody products 研发新抗体产品时的分析方法


Feature 特点	Analysis method 分析方法
Ig class and subclass 抗体类型及亚型	Commercial lateral flow rapid test 工业用快速免疫层析测试
Stability 稳定性	Accelerated stability tests at elevated temperatures 热加速稳定性测试
Immunoreactivity 免疫学活性	Fluoroimmunoassay, direct coating and/or sandwich assay 免疫荧光测试，直接包被和/或夹心测试
Pair testing - with purified antigen and clinical samples 配对测试（纯化抗原、患者样本）	Sandwich-FIA, lateral flow when applicable 根据可应用性，夹心法免疫荧光测试，免疫层析
Purity 纯度	CE (Capillary electrophoresis) 毛细管电泳
Homogeneity 同质性	HPLC-SEC (size exclusion chromatography) HPLC-SEC(分子筛色谱)
Charge isoform profile 电荷亚型分布	cIEF (capillary isoelectric focusing) (毛细管等电点聚焦)
Kinetic constants 动力学常数	Surface Plasmon Resonance (SPR) 表面等离子体共振 (SPR)
Concentration 浓度	Spectrophotometer (Abs 280 nm) 分光光度计 (280nm)

Confidential - distribution prohibited without permission.

Example: Antibodies are developed to meet IVD industry requirements 范例：根据IVD企业需求研发抗体

Patient samples were measured with a sandwich fluoroimmunoassay using Medix Biochemica's anti-NT-proBNP mouse monoclonal antibodies 1306 (capture) and a newly developed 1309 (label). The results showed excellent correlation with a reference NT-proBNP method:

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Confidential - distribution prohibited without permission.

Validated antibody production methods

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
- Optimized, serum and protein free, industrial-scale *in vitro* production methods
经优化的无血清和蛋白的工业级体外生产方法
- Hollow fiber bioreactors or Wave tanks
中空纤维生物反应器或波浪摇床
- Antibody purification by Protein A affinity chromatography
采用蛋白A亲和色谱进行抗体纯化

Confidential - distribution prohibited without permission.

Hollow fiber cell culture 中空纤维细胞培养

- Best suited for production scale 1 g →
最适生产规模 1g
 - Perfusion cultivation, constant fresh media delivery 充填式培养，新鲜培养基的持续输送
 - Advantage: high antibody concentration in the supernatant, long protein production times
优势：上清液中存在高浓度抗体，可长时间持续生产
 - Linearly scalable; bioreactors available in different sizes 可线性扩展，不同大小的生物反应器可供选择

Confidential - distribution prohibited without permission.

Validated antibody quality control methods

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
QC information can be found on Product Specification Sheets and on Certificates of Analysis. Product Specifications are available online: www.medixbiochemica.com

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Confidential - distribution prohibited without permission.

Immunoreactivity test 免疫学活性检测

- Each new antibody lot is compared against a reference lot using a wide antibody concentration range in a fluoroimmunoassay
通过免疫荧光法，每个新批号抗体都需要在很宽的抗体浓度范围与参考批号进行比较。
- The antibody concentration showing 50% of the maximum binding is calculated
计算50%最大结合度的抗体浓度
- Sample / reference result must within 80-120%
样品/参考抗体的结果必须在80-120%之间
- The tolerance limits are dictated by the FIA method
允许范围由免疫荧光方法学决定。


Confidential - distribution prohibited without permission.

Medix Biochemica

Example: High batch-to-batch consistency 范例：高度的批间一致性

Immunoreactivity quality control results for MedixMAB SAA 2201,
51 batches manufactured and analyzed in December 2011 - February 2018
MedixMAB SAA 2201的免疫学活性质检结果，从2011年12月至2017年2月的51个生产批次

Confidential - distribution prohibited without permission.

New products in the pipeline for 2018 2018年计划新产品

- **Monoclonal antibodies** against
 - NT-proBNP – new clones to pair with mAb 1306
 - Cystatin C
 - Calcitonin
 - T3
 - Cardiac Troponin T
 - *Helicobacter pylori*
 - *Clostridium difficile* GDH
 - Rotavirus
 - ST2
 - Galectin-3
- **Antigens**
 - Influenza A and B nucleoproteins
 - Galectin-3

Confidential - distribution prohibited without permission.

Confidential - distribution prohibited without permission.

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Lluis M. Martinez, PhD
Chief Scientific Officer

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