


25(OH) Vitamin D2/D3

(Mouse monoclonal Antibody anti Vitamin D2 / D3)



PRODUCT ANALYSIS SHEET

Cat 53.197.16

Clone:	LMBP7012CB
Batch:	B9I2P1
Size	2 mg / 290 µl
Format:	Purified , Liquid (Phosphate Buffer)
Concentration:	6.91 mg/ml in PBS 0.1M pH 7.4
Preservative:	NaN3 (0.05%)
Purification:	Purified by Protein A affinity chromatography
Isotype:	Ig G1 kappa
Immunogen:	Vitamin D derivate - BSA
Mouse / fusion partners:	Balb/c (Charles River)/ NSO
Specificity :	Recognizes natural Vitamin D2 and Vitamin D3
Application:	Elisa
Storage & Stability:	2-8 /10 years
Update:	

Date & Signature

DIAsource ImmunoAssays S.A.

Tel +32 10 84 99 11 | Fax +32 10 84 99 90 | www.diasource.be | Technical Inquiries: tech.support@diasource.be
To order: Belgium Free Phone 0800 159 59 | BeLux Tel +32 (0) 10 84 99 00 | Fax +32 (0) 10 84 99 96 | ordering@diasource.be
France Free Phone 0800908 443 | Free Fax 0800 902 588 | ordering@diasource.be