

EastCoast Bio, Inc.

211 Wells Street | P.O. Box 489 | North Berwick, ME 03906

P 207.676.7639 | F 207.676.7658 | E info@eastcoastbio.com

Product Data: SeaBlock w/ Tris

Catalog #	PP81-Tris
Lot #	Y2632
Description	A non-mammalian blocking reagent supplied in TBS for diagnostic systems involving membrane or other solid phase material. This new reagent imparts superior stability and flow characteristics, while utilizing proteins that do not interact with mammalian antibodies.
Format	Liquid
pH	7.6 +/- 0.2
Expires	3/15/2020
Preservative	0.09% Sodium Azide
Dilution/Recommended Use	Ready to use formulation. Use at 1x concentration as well as further dilutions which are dependent upon specific user's application.
Applications	ELISA, WB
Storage	2-8°C. Do not freeze.
Note	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.