

Product Data: SeaBlock w/ PBS

Catalog #	PP81
Lot #	Y0712
Description	A non-mammalian blocking reagent supplied in PBS for diagnostic systems involving membrane or other solid phase material. This new reagent is optimized for superior performance in super sensitive assay systems, while utilizing proteins that do not interact with mammalian antibodies.
Format	Liquid
pH	7.4 +/- 0.1
Expires	3/1/2020
Preservative	0.09% Sodium Azide
Dilution/Recommended Use	Ready to use formulation. Use at 1x concentration as well as further dilutions which are dependent upon specific user's application.
Applications	ELISA, WB
Storage	2-8°C. Do not freeze.
Note	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.