

## Product Data: SeaBlock Serum Free w/ PBS

<b>Catalog #</b>	PP83A
<b>Lot #</b>	Z0744
<b>Description</b>	A non-mammalian blocking reagent supplied in PBS for diagnostic systems involving membrane or other solid phase material. This new reagent imparts superior stability and flow characteristics, while utilizing proteins that do not interact with mammalian antibodies.
<b>Format</b>	Liquid
<b>pH</b>	7.4 +/- 0.2
<b>Expires</b>	9/1/2020
<b>Preservative</b>	0.09% Sodium Azide
<b>Dilution/Recommended Use</b>	Ready to use formulation. Use at 1x concentration as well as further dilutions which are dependent upon specific user's application.
<b>Applications</b>	Nitrocellulose in Lateral Flow
<b>Storage</b>	2-8°C. Do not freeze.
<b>Note</b>	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.