

EastCoast Bio, Inc.

211 Wells Street | P.O. Box 489 | North Berwick, ME 03906

P 207.676.7639 | F 207.676.7658 | E info@eastcoastbio.com

Product Data: Balance Blocking Buffer w/ PBS

Catalog #	PP91
Lot #	Z0933
Description	A universal blocking reagent supplied in PBS for diagnostic and research systems. This new reagent contains a unique blend of non-mammalian proteins providing superior blocking in super sensitive assay systems.
Format	Liquid
pH	7.4 +/- 0.1
Preservative	0.09% Sodium Azide
Dilution/Recommended Use	Ready to use formulation. Use at 1x concentration as well as further dilutions which are dependent upon specific user's application.
Applications	ELISA, WB, Lateral Flow
Storage	2-8°C. Do not freeze.
Note	This product is intended for research and manufacturing uses only. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.