

Immunology Consultants Laboratory, Inc.

15862 SW 72nd Avenue, #150
Portland, OR 97224
T 503-747-2454
F 503-747-2544

www.icllab.com



ANTIBODIES ARE OUR BUSINESS

PRODUCT PROFILE

PRODUCT: Goat anti-Human Fibronectin - Antiserum

CATALOG NUMBER: GFCT-80

LOT NUMBER: 23

QUANTITY: 16,000.0mL (16 bottles @ 1,000.0mL/Bottle)

PHYSICAL STATE: Liquid

PRESERVATIVE: 0.1% Sodium Azide

STORAGE: 4° C

EXPIRATION: 24 months from shipment date

PROCEDURE: Goats were immunized with purified human fibronectin. Serum was collected and delipidated, then solid phase absorbed to ensure specificity.

SPECIFICATIONS: This antiserum shows a precipitin band by immunoelectrophoresis against purified human fibronectin, human plasma and serum.

Total protein value as measured by refractometer = 71 mg/mL.

Vialed Reagents
For In Vitro Laboratory Use Only!!
Not for Human or Animal Consumption!!
Store at 4C

FOR RESEARCH AND MANUFACTURING ONLY; NOT FOR USE IN DIAGNOSTIC APPLICATIONS.
FOR IN VITRO USE ONLY, NOT FOR HUMAN OR ANIMAL CONSUMPTION

CERTIFICATE OF ANALYSIS VERIFICATION:
PRODUCT OF THE UNITED STATES OF AMERICA

DATE: 07/19/19

GFCT80